



Order online: concord-surgical.com
 Call toll free 1-877-330-3335
 Fax orders 1-888-897-8081

OCT FLYER

SALE



Savings on surgical instruments and much more

Concord

Oral Medication Kits



Everything your practice needs for efficient oral med dispensing.
Price is for bag of 12 kits

Each kit contains :

- One 1 oz round blue bottle
- One child resistant closure
- One push-in bottle adapter cap
- One 1oz wrapped oral syringe

Code: OEL-OMK

Reg \$22

SALE \$19.99



Concord

Locking Cancellous Screws - 4.0mm

Exclusive to Concord!. Our 4.0mm locking cancellous screws have been fitted with the 3.5mm locking heads for cortical screws, making the head compatible with 3.5mm TPLO and DCP locking plates. Lowers the possibility of "loosening", screws remains locked in place. Lowers possibility of post-op failures



4.0mm shaft

316L implant stainless steel

3.5mm locking head

Self tapping, double fluted

Hex-head, fully threaded

Reg \$15

SALE

\$12.99

Call or goto concord-surgical.com for codes

20% OFF



Concord

Lister Bandage Scissors

German stainless steel 

Size	Code	Reg price	Sale price
9cm	OEL-27-100	\$28	\$22.40
11cm	Oel-27-102	\$30	\$24.00
14cm	OEL-27-104	\$34	\$27.20

Concord

Adson Tissue Forceps 1:2 teeth

German stainless steel 

Reg \$20

SALE \$15.99



Code: OEL-15-330

Rochester Carmalt Haemostatic Forceps

German stainless steel 

Reg \$40

Code: OEL-17-410 - STR 16cm

Code: OEL-17-420 - CRV 16cm

SALE \$32



Concord

Snooks Spay Hook

German stainless steel 

Reg \$20

SALE \$15.99



Code: OJO-J112A

20% OFF

Mayo Scissors German stainless steel



Designed to cut heavy fascia and sutures

20% OFF

Concord

Size /description	Code	Reg price	Sale price
14cm straight	OEL-13-180	\$28	\$22.40
14cm curved	OEL-13-190	\$28	\$22.40
17cm straight	OEL-13- 184	\$30	\$24.00
17cm curved`	OEL-13-194	\$30	\$24.00
17cm straight t-carbide	OEL-13-184T	\$115	\$92.00
17cm curved t- carbide	OEL-13-194T	\$120	\$96.00

Gelpi Retractors

15% OFF

German stainless steel

size	code	reg price	sale price
9cm	OEL-21-461	100	\$85
11cm	OEL-21-462	\$105	\$89
13cm	OEL-21-460	\$113	\$96



Concord



Disposable Surgical Drape

Reg \$150

SALE \$135

Blue, soft to the touch, latex free material
Autoclavable: 100yd roll x 38inches wide

Code: DSD100



Winged Elevator Sets Coloured, Stubby Handles



French stainless steel

Sets of 4 or 5 with choice of case, leather or autoclavable stainless steel

20% OFF

Concord

description	Code	Reg price	Sale price
Set 4 leather case	MED-733SET4	\$169	\$135
Set 4 stainless case	MED-733SET4SSC	\$249	\$199
Set 5 leather case	MED-733SET5	\$199	\$159
Set 5 stainless case	MED-733SET5SSC	\$289	\$229



Riester

Fortelux[®] Diagnostic Penlight



Made in Germany

2.2V vacuum lamp

Disinfectable aluminum

Professional medical quality

Reg \$30

SALE \$19.99

Bishop Harmon Forceps

Light action forceps ideal for grasping of delicate tissues in eye surgery

German stainless steel

15% OFF



Concord

Size /description	Code	Reg price	Sale price
9cm 1:2 teeth	OEL-37-642	\$40	\$34
9cm -serrated tip	OEL-37-640	\$40	\$34

Cold Pack



Reg \$192

Concord
INTEGRA
Miltex

SALE

Miltex germicide tray
Olsen Hegar needle holder
Scalpel handle #3
Curved mosquito forceps
Mayo scissors- straight
Adson 1:2 thumb forceps

SALE \$169.99

To be removed from email or fax list, please email concordvet supply@telus.net or call our toll free line 1-877-330-3335